


# JOBST® FarrowWrap®

## A-G Thigh High

Ready-to-wear order form and prescription request

Customer Services on: **031 710 8111**  
 Email order to Customer Services on: [medical.za@essity.com](mailto:medical.za@essity.com)

Date: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Patient Name: \_\_\_\_\_ DoB: \_\_\_\_\_  
 Measured By: \_\_\_\_\_ Tel: \_\_\_\_\_ Email: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_ Invoice Address: \_\_\_\_\_



## Measuring instructions

### Thighpiece - 3 simple measuring points

1. Measure 5cm above the patella **E** and the groin **G** circumference (as shown on diagram)
2. Use chart below to determine size at the top of the column (if 2 measurements do not fit in same column the smaller size should accommodate reducing oedema)
3. Measure the length of the thigh medially between **E** and **G** following the contour of the leg. Use the sizing chart to select short, regular or tall to determine size

### Legpiece - 3 simple measuring points

1. Measure the widest calf **C** and narrowest ankle **B** circumference (as shown on diagram)
2. Use chart below to determine size at the top of the column (if 2 measurements do not fit in same column the smaller size should accommodate reducing the oedema)
3. Measure the length of the leg posteriorly from the floor to **D** measurement (2cm below the crease of the knee). Follow the contour of the leg. Use the sizing chart to select regular or tall to determine size

### Footpiece - 2 simple measuring points

1. Measure **A'** circumference (as shown on diagram)
2. Use chart below to determine size at the top of the column
3. Measure the **X** length of the foot. Use the sizing chart to select regular or long to determine size

	X Small	Small	Medium	Large	X Large
<b>Thighpiece (including kneepiece)</b>					
<b>G</b> Groin	60-70cm	65-75cm	70-80cm	75-90cm	85-100cm
<b>E</b> Above knee	47-52cm	53-59cm	60-66cm	67-73cm	74-80cm
<b>E-G</b> Short	16-20cm	16-20cm	16-20cm	16-20cm	16-20cm
<b>E-G</b> Regular	21-26cm	21-26cm	21-26cm	21-26cm	21-26cm
<b>E-G</b> Tall	27-31cm	27-31cm	27-31cm	27-31cm	27-31cm
NOTE: for <b>E-G</b> measure medially and follow contours					
<b>Legpiece</b>					
<b>C</b> Mid-Calf	36-43cm	42-50cm	48-58cm	53-63cm	58-68cm
<b>B</b> Ankle	21-25cm	25-30cm	30-36cm	36-42cm	42-50cm
<b>A-t</b> Regular	35-39cm	37-41cm	39-43cm	41-45cm	42-46cm
<b>A-t</b> Tall	40-44cm	42-46cm	44-48cm	46-50cm	47-51cm
NOTE: for <b>A-D</b> measure posteriorly and follow contours					
<b>Footpiece</b>					
<b>A'</b> Mid-foot	22-24cm	25-27cm	28-30cm	31-34cm	35-40cm
<b>X</b> Regular	15-17cm	16-19cm	18-20cm	19-22cm	21-24cm
<b>X</b> Tall	18-20cm	20-22cm	21-23cm	23-25cm	25-27cm
NOTE: for <b>X</b> measure from the heel to the base of the little toe					

**Thighpiece**

**Size:**  X Small  Small  Medium  Large  X Large

**Length E-G:**  Short  Regular  Tall **Fabric:**  Lite  Classic  Strong

PIP Code:  -

Code:  -

Quantity Required:

Repeat prescription required every  months

**Legpiece**

**Size:**  X Small  Small  Medium  Large  X Large

**Length A-D:**  Regular  Tall **Fabric:**  Lite  Classic  Strong

PIP Code:  -

Code:  -

Quantity Required:

Repeat prescription required every  months

**Footpiece**

**Size:**  X Small  Small  Medium  Large  X Large

**Length X:**  Regular  Long **Fabric:**  Lite  Classic  Strong

PIP Code:  -

Code:  -

Quantity Required:

Repeat prescription required every  months

**JOBST® FarrowHybrid®**

**Size:**  Small  Medium  Large  X Large

**Width:**  Regular\*  Wide\*\* **Compression:** 20-30mmHg -  Taupe  Black\*\*\*

PIP Code:  -

Code:  -

Quantity Required:

Repeat prescription required every  months

A pair of appropriately sized liners will be provided free of charge with each JOBST® FarrowWrap® lower limb piece.  
 Below knee comes with silver liners, thigh high comes with TG Soft liners.

Customer Service: **031 710 8111**  
 Email: [medical.za@essity.com](mailto:medical.za@essity.com) Website: [www.jobst.co.za](http://www.jobst.co.za)

Note: 3cm allowance on either side of measurement.

\*Regular: 20-27 cm **A'** \*\*Wide: 28-35 cm **A'** \*\*\*Available in wide variant only

# JOBST® FarrowWrap® Order Codes

## JOBST® FarrowWrap® Lite (20-30mmHg)

Code	Size	Length	Colour	PIP Code
<b>Thighpiece (including kneepiece)</b>				
76655-05	Extra Small	Short	Beige	360-6084
76655-00	Extra Small	Regular	Beige	360-6134
76655-10	Extra Small	Tall	Beige	360-6183
76655-06	Small	Short	Beige	360-6092
76655-01	Small	Regular	Beige	360-6142
76655-11	Small	Tall	Beige	360-6191
76655-07	Medium	Short	Beige	360-6100
76655-02	Medium	Regular	Beige	360-6159
76655-12	Medium	Tall	Beige	360-6209
76655-08	Large	Short	Beige	360-6118
76655-03	Large	Regular	Beige	360-6167
76655-13	Large	Tall	Beige	360-6217
76655-09	Extra Large	Short	Beige	360-6126
76655-04	Extra Large	Regular	Beige	360-6175
76655-14	Extra Large	Tall	Beige	360-6225
<b>Legpiece</b>				
76654-00	Extra Small	Regular	Beige	360-6233
76654-05	Extra Small	Tall	Beige	360-6282
76654-01	Small	Regular	Beige	360-6241
76654-06	Small	Tall	Beige	360-6290
76654-02	Medium	Regular	Beige	360-6258
76654-07	Medium	Tall	Beige	360-6308
76654-03	Large	Regular	Beige	360-6266
76654-08	Large	Tall	Beige	360-6316
76654-04	Extra Large	Regular	Beige	360-6274
76654-09	Extra Large	Tall	Beige	360-6324
<b>Footpiece</b>				
76652-05	Extra Small	Regular	Beige	360-6332
76652-00	Extra Small	Long	Beige	360-6381
76652-06	Small	Regular	Beige	360-6340
76652-01	Small	Long	Beige	360-6399
76652-07	Medium	Regular	Beige	360-6357
76652-02	Medium	Long	Beige	360-6407
76652-08	Large	Regular	Beige	360-6365
76652-03	Large	Long	Beige	360-6415
76652-09	Extra Large	Regular	Beige	360-6373
76652-04	Extra Large	Long	Beige	360-6423



## JOBST® FarrowWrap® Classic (30-40mmHg)

Code	Size	Length	Colour	PIP Code
<b>Thighpiece (including kneepiece)</b>				
76658-06	Extra Small	Short	Beige	360-6431
76658-00	Extra Small	Regular	Beige	360-6480
76658-12	Extra Small	Tall	Beige	360-6530
76658-07	Small	Short	Beige	360-6449
76658-01	Small	Regular	Beige	360-6498
76658-13	Small	Tall	Beige	360-6548
76658-08	Medium	Short	Beige	360-6456
76658-02	Medium	Regular	Beige	360-6506
76658-14	Medium	Tall	Beige	360-6555
76658-09	Large	Short	Beige	360-6464
76658-03	Large	Regular	Beige	360-6514
76658-15	Large	Tall	Beige	360-6563
76658-10	Extra Large	Short	Beige	360-6472
76658-04	Extra Large	Regular	Beige	360-6522
76658-16	Extra Large	Tall	Beige	360-6571
<b>Legpiece</b>				
76657-00	Extra Small	Regular	Beige	360-6589
76657-05	Extra Small	Tall	Beige	360-6639
76657-01	Small	Regular	Beige	360-6597
76657-06	Small	Tall	Beige	360-6647
76657-02	Medium	Regular	Beige	360-6605
76657-07	Medium	Tall	Beige	360-6654
76657-03	Large	Regular	Beige	360-6613
76657-08	Large	Tall	Beige	360-6662
76657-04	Extra Large	Regular	Beige	360-6621
76657-09	Extra Large	Tall	Beige	360-6670
<b>Footpiece</b>				
76656-05	Extra Small	Regular	Beige	360-6688
76656-00	Extra Small	Long	Beige	360-6738
76656-06	Small	Regular	Beige	360-6696
76656-01	Small	Long	Beige	360-6746
76656-07	Medium	Regular	Beige	360-6704
76656-02	Medium	Long	Beige	360-6753
76656-08	Large	Regular	Beige	360-6712
76656-03	Large	Long	Beige	360-6761
76656-09	Extra Large	Regular	Beige	360-6720
76656-04	Extra Large	Long	Beige	360-6779



## JOBST® FarrowWrap® Strong (30-40mmHg)

Code	Size	Length	Colour	PIP Code
<b>Thighpiece (including kneepiece)</b>				
76649-05	Extra Small	Short	Beige	382-0222
76649-00	Extra Small	Regular	Beige	382-0172
76649-10	Extra Small	Tall	Beige	382-0123
76649-06	Small	Short	Beige	382-0255
76649-01	Small	Regular	Beige	382-0206
76649-11	Small	Tall	Beige	382-0156
76649-07	Medium	Short	Beige	382-0248
76649-02	Medium	Regular	Beige	382-0198
76649-12	Medium	Tall	Beige	382-0149
76649-08	Large	Short	Beige	382-0230
76649-03	Large	Regular	Beige	382-0180
76649-13	Large	Tall	Beige	382-0131
76649-09	Extra Large	Short	Beige	382-0214
76649-04	Extra Large	Regular	Beige	382-0164
76649-14	Extra Large	Tall	Beige	382-0115
<b>Legpiece</b>				
76648-00	Extra Small	Regular	Beige	382-0073
76648-05	Extra Small	Tall	Beige	382-0024
76648-01	Small	Regular	Beige	382-0107
76648-06	Small	Tall	Beige	382-0057
76648-02	Medium	Regular	Beige	382-0099
76648-07	Medium	Tall	Beige	382-0040
76648-03	Large	Regular	Beige	382-0081
76648-08	Large	Tall	Beige	382-0032
76648-04	Extra Large	Regular	Beige	382-0065
76648-09	Extra Large	Tall	Beige	382-0016
<b>Footpiece</b>				
76646-04	Extra Small	Regular	Beige	381-9976
76646-09	Extra Small	Long	Beige	381-9927
76646-05	Small	Regular	Beige	382-0008
76646-00	Small	Long	Beige	381-9950
76646-06	Medium	Regular	Beige	381-9992
76646-01	Medium	Long	Beige	381-9943
76646-07	Large	Regular	Beige	381-9984
76646-02	Large	Long	Beige	381-9935
76646-08	Extra Large	Regular	Beige	381-9968
76646-03	Extra Large	Long	Beige	381-9919

## JOBST® FarrowHybrid®\* - supplying compression to the foot and ankle

Code	Size	Colour	PIP Code
<b>ADI (20-30mmHg)</b>			
76664-00	Small (shoe size 2-5)	Beige	408-1626
76664-01	Medium (shoe size 6-8)	Beige	408-1634
76664-02	Large (shoe size 9-11)	Beige	408-1642
76664-03	Extra Large (shoe size 12-13)	Beige	408-1659
76664-12	Small - wide (shoe size 2-5)	Beige	408-1667
76664-13	Medium - wide (shoe size 6-8)	Beige	408-1675
76664-14	Large - wide (shoe size 9-11)	Beige	408-2384
76664-15	Extra Large - wide (shoe size 12-13)	Beige	408-1683
76664-16	Small - wide (shoe size 2-5)	Black	408-6765
76664-17	Medium - wide (shoe size 6-8)	Black	408-6757
76664-18	Large - wide (shoe size 9-11)	Black	408-6740
76664-19	Extra Large - wide (shoe size 12-13)	Black	408-6732



\*Do not apply JOBST® FarrowHybrid® if calf circumference exceeds 60cm. Consider a JOBST® FarrowWrap® footpiece in this circumstance.

**Prescription request**  
If requesting a prescription please complete order form and return to the prescriber. The codes for specific products can be found on this form

By completing this order form, you are confirming that you are aware of your obligation to obtain consent from the patient on the processing of their data for the production of their JOBST® compression garment in accordance with the General Data Protection Regulation. For more information and our Privacy Policy, visit [www.jobst.co.za](http://www.jobst.co.za)