

READY-TO-WEAR ORDER FORM LOWER EXTREMITY CUSTOMER SERVICE: 031 710 8111

Date: _____ Purchase Order No.: _____

Patient Name: _____

Account number and name: _____

Fax number: _____

Circumference and length measurements in centimeters (cms)

NOTE:
For A-D and D-E, measure posteriorly and follow leg contour.
For E-G measure medially and follow contours.

E-G

Left:

Right:

D-E
Note: only for CM standard knee


Left:

Right:

A-D

Left:

Right:



Remarks

STEPS: 1. Measure the patient at all ■ measuring points. 2. Determine if patient fits within ready-to-wear (RTW) sizing range by using the size chart. 3. If patient does not fit within a RTW, take the additional measurements to order a custom-made (CM) piece and send to Essity. 4. It is possible to mix and match RTW and CM pieces for a patient.

*To reduce the risk of error, orders will only be placed if the product code is indicated.

	XSmall	Small	Medium	Large	XLarge
Legpiece (except JOBST® FarrowWrap® 4000)					
C Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-68 cm
B Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm
A-D Regular	33-37 cm	35-39 cm	37-41 cm	39-43 cm	39-43 cm
A-D Tall	38-41 cm	40-43 cm	42-45 cm	44-47 cm	44-47 cm
Colour <input checked="" type="checkbox"/> Tan			Side <input type="checkbox"/> Left <input type="checkbox"/> Right		
Product code*	<input type="text"/>				
Thighpiece					
G Groin	60-70 cm	67-75 cm	70-80 cm	75-90 cm	85-100 cm
E Above patella	47-52 cm	53-59 cm	60-66 cm	67-73 cm	74-80 cm
Colour <input checked="" type="checkbox"/> Tan			Side <input type="checkbox"/> Left <input type="checkbox"/> Right		
Product code*	<input type="text"/>				
Footpiece					
A' Mid-foot	22-24 cm	25-27 cm	28-30 cm	31-34 cm	35-40 cm
X Regular	16-17 cm	17-18 cm	19-20 cm	20-21 cm	22-23 cm
X Tall	18-19 cm	19-20 cm	21-22 cm	22-23 cm	24-25 cm
Colour <input checked="" type="checkbox"/> Tan			Side <input type="checkbox"/> Left <input type="checkbox"/> Right		
Product code*	<input type="text"/>				
Liners*					
C Silver (standard)		19-48 cm		48-70 cm	
C TG® Soft		<40 cm	40-70 cm	70-120 cm	
FarrowHybrid™ Shoe size Men		3-6	7-9	10-12	13-15
Foot Compression Shoe size Women		4-7	8-10	11-13	14-16
FarrowHybrid™ ADI (20-30 mmHg) <input type="checkbox"/>			FarrowHybrid™ ADII (30-40 mmHg) <input type="checkbox"/>		
*Supplied in pairs					
Notes for Legpiece: Please use JOBST® FarrowWrap® 4000 order form.					
Notes for FarrowHybrid™: Calf circumference < 60 cm. Standard foot A' = 20-27cm Wide foot A' = 28-35cm					